



MiniMed Flex™ system

How to prescribe



Nothing comes close to freedom you can feel with fewer changes[†]

The MiniMed Flex™ system

- Meal Detection™ technology[§]
- Flexible targets as low as 100 mg/dL
- Recommended settings: 100 mg/dL target with 2-hr AIT
- Missed bolus forgiveness[§]
- Real-world results[¶] using SmartGuard™ feature showed 80% Time in Range¹ with recommended settings^{¶¶}

Life with Flexibility

- Small and discreet
- Full smartphone-controlled operation
- Clean and modern waterproof^{¶¶} design
- 300U reservoir
- Up to 7-day infusion set
- 30-minute charge provides up to 7 days of use

Learn more about the MiniMed Flex system:



Download the MiniMed Simulator App
Experience the ease of onboarding – use the app as a safe space for new users to practice workflows.



Listen to the Podcast
Subscribe to listen to the latest Advancing the dose podcast.



Visit professional.minimed.com for the latest clinical resources, to request an in-service, and to register for upcoming webinars and livestreams



Help your patients get started by taking the next step

How to Prescribe?



Submit your prescription to MiniMed through Parachute Health.



Submit your prescription to MiniMed via fax:

Fax: 818-576-6641

If we need any additional information, we will contact your office.



Scan to download CMN.

Process to Prescribe

Directions section of the ePrescribe screen or on written prescriptions:

MiniMed Flex™ pump

Include the following:

- Indicate SmartGuard™ Target: 100 (recommended as clinically appropriate),^Δ 110, and 120 mg/dL
Note: System will default to 100 mg/dL.
- Indicate Active Insulin Time: 2 hours (recommended)
- SmartGuard™ Warmup/Manual Mode: Suspend before low (recommended) or Suspend on low

Extended™ infusion set

Include the following:

- Change infusion set & reservoir every 7 days (recommended)

Extended™ Reservoir

Include the following: Total Daily Dose

Simplera Sync™ sensor

- Indicate change sensor: Up to 6 days, followed by a 24 hour grace period (recommended)

How can your patients order the MiniMed Flex™ system?

Existing Patients:



Encourage your patients to check their upgrade eligibility on [Diabetes.shop](https://diabetes.shop).

New Patients:



Visit minimed.com or reach out to your local territory team.



[Diabetes.shop](https://diabetes.shop) is where they can complete required training.

† Based on usage of 4 Simplera Sync sensors for 6 days each and 4 Extended Infusion Sets for 7 days each in a one-month period. More insertions may be required if Simplera Sync sensor or Extended infusion set requires replacement before specified wear time. Refers to auto correct, which provides bolus assistance. Can deliver all correction doses automatically without user interaction, feature can be turned on and off.

‡ Refers to SmartGuard™ feature. Individual results may vary.

§ Taking a bolus 15 – 20 min before a meal provides significant improvement in post meal control. We recommend taking a meal bolus.

¶ Due to inherent real-world study limitations, caution is advised when attempting to extrapolate these results to new patients. There could be significant differences.

At the time of manufacture and when the reservoir and tubing are properly inserted, the pump is waterproof. It is protected against the effects of being underwater to a depth of up to 8 feet (2.4 meters) for 1 hour and 45 minutes if you choose to stay connected while in water. This classified as IPX8 rating. See user guide for more details. The sensor is water-resistant at 8 feet (2.4 meters) for up to 30 minutes. CGM readings may not be transmitted from the CGM to the pump while in water.

Δ Glucose targets should be set following a consultation with your healthcare provider.

◊◊ The SmartGuard™ feature in the MiniMed Flex is equivalent to the MiniMed 780G system.

1. Arrieta A, et al. Diabetes Obes Metab. 2022;10.1111/dom.14714

Important Safety information for MiniMed Flex™ system:

MiniMed Flex™ System is for type 1 ages 7 and over and type 2 ages 18 and over. Prescription required. **WARNING: Do not use SmartGuard™ feature for people who require less than 8 units or more than 250 units of insulin/day.** For additional safety information and user guides, see <https://bit.ly/MiniMedRisks>.

US-DBA-2600624 ©2026 Medtronic MiniMed, Inc. MiniMed and MiniMed logo are trademarks of Medtronic MiniMed, Inc.

™™Third-party brands are trademarks of their respective owners.